














Biotrak®	Biotrak Mini  16.0 mm	Biotrak Standard  16.0 mm	Biotrak Pin  20.0 mm	Biotrak Helical Nail  30.0 mm
<b>Biotrak Part No.</b>	301700XX-S	301700XX-S	301703XX-S	55-00XX-S
<b>Lengths</b>	16.0 mm–24.0 mm (2.0 mm increments)	16.0 mm–24.0 mm (2.0 mm increments)	20 mm, 30 mm, 40 mm	20 mm, 30 mm, 40 mm
<b>Tip Diameters</b>	3.2 mm tip	3.6 mm tip	2.0 mm tip	2.5 mm tip
<b>Tail Diameters</b>	3.5 mm–3.7 mm tail	4.3 mm–4.7 mm tail	3.0 mm tail	3.2 mm tail
<b>Guide Wire</b>	80100101	80100100	80100102 (Single Trocar) 80100103 (Double Trocar)	.035X8 DT 80-0470 .035X8 ST 80-0471 .045X8 DT 80-0328 .045X8 ST 80-0329
<b>Guide Wire Diameter</b>	.035" (.9 mm)	.045" (1.1 mm)	.035" (.9 mm)	.035" (.9 mm) .045" (1.1 mm)
<b>Drill Part No.</b>	30100151	30100101	30100351	80-0327 (2.5mm) 80-0469 (2.0mm)
<b>Driver Part No.</b>	30100153	30100103	30100353 (Plunger)	80-0323 (Plunger)


AcuTwist®	AcuTwist  10.0 mm
<b>AcuTwist Part No.</b>	AI-00XX-S (non-cannulated)
<b>Lengths</b>	10.0 mm–30.0 mm (2.0 mm increments)
<b>Tip Diameters</b>	1.5 mm tip
<b>Tail Diameters</b>	2.0 mm tail
<b>Guide Wire</b>	WS-1106ST
<b>Guide Wire Diameter</b>	.045" (1.1 mm)
<b>Tap Part No.</b>	AI-NG30



Acutrak 2®	Acutrak 2: Micro  8.0 mm	Acutrak 2: Mini  16.0 mm	Acutrak 2: Standard  16.0 mm	Acutrak 2: 4.7  20.0 mm	Acutrak 2: 5.5  25.0 mm	Acutrak 2: 7.5  40.0 mm
<b>Acutrak 2 Part No.</b>	AT2-CXX-S and AT2-CXX	AT2-MXX-S and AT2-MXX	AT2-SXX-S and AT2-SXX	30-0XXX and 30-0XXX-S	30-00XX and 30-00XX-S	30-0XXX and 30-0XXX-S
<b>Lengths</b>	8.0 mm–14.0 mm (1.0 mm increments), 16.0 mm–30.0 mm (2.0 mm increments)	16.0 mm–30.0 mm (2.0 mm increments)	16.0 mm–34.0 mm (2.0 mm increments)	20.0 mm–30.0 mm (2.0 mm increments) 30.0 mm–50.0 mm (5.0 mm increments)	25.0 mm–60.0 mm (5.0 mm increments)	40.0 mm–120.0 mm (5.0 mm increments)
<b>Tip Diameters</b>	2.5 mm tip	3.5 mm tip	4.0 mm tip	4.5 mm tip	5.2 mm tip	7.0 mm tip
<b>Tail Diameters</b>	2.8 mm tail	3.6 mm tail	4.1 mm tail	4.7 mm tail	5.5 mm tail	7.5 mm tail
<b>Guide Wire</b>	WS-0906ST, 80-1524 Single Trocar, 80-1525 Double Trocar	WS-1106ST	WS-1407ST	80-0950	80-0950 80-0413 in AT2-5.5 System	80-0970 non-threaded 80-0971 threaded
<b>Guide Wire Diameter</b>	.035" (.9 mm)	.045" (1.1 mm)	.054" (1.4 mm)	.062" (1.6 mm)	.062" (1.6 mm)	.094" (2.4 mm)
<b>Hex Driver</b>	HT-0915	HT-1120	HT-1725	80-0958 (cannulated) 80-0959 (solid)	80-0958 (cannulated) 80-0959 (solid)	80-0978 (cannulated) 80-0979 (solid)
<b>Hex Size</b>	1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.0 mm	4.0 mm
<b>Profile Drill Part No.</b>	AT2-1509	AT2M-1813	AT2-2515	80-0945	80-0955, 80-0055 in AT2-5.5 System	80-0975
<b>Profile Drill Tip Diameter</b>	1.8 mm tip	2.5 mm tip	3.1 mm tip	2.8 mm tip	4.0 mm tip	4.5 mm tip
<b>Profile Drill Tail Diameter</b>	2.9 mm tail	3.4 mm tail	4.1 mm tail	4.7 mm tail	5.5 mm tail	6.9 mm tail
<b>Profile Drill Length</b>	4.6 mm	8.2 mm	8.2 mm	16.2 mm (to 2mm band)	14.5 mm (to 2 mm band)	24.3 mm (to 2mm band)
<b>Long Drill Part No.</b>	80-0100, 80-1522	AT2M-L1813	AT2-L2515	80-0946	80-0956, 80-0056 in AT2-5.5 System	80-0976
<b>Long Drill Diameter</b>	1.9 mm	2.7 mm	3.2 mm	3.3 mm	4.0 mm	5.3 mm
<b>Long Drill Length</b>	20.0 mm, 36.0 mm	30.0 mm	34.0 mm	50.0 mm	60.0 mm	120.0 mm



# Screw & Pin Fixation Reference Chart

<b>Acutrak® Fusion</b>	<b>Acutrak Fusion</b>  <b>14.0 mm</b>
<b>Acutrak Part No.</b>	ATF-XXX-S-30-0033-S (non-cannulated)
<b>Lengths</b>	14.0 mm–24.0 mm (2.0 mm increments) 27.0 mm, 30.0 mm, 32.0 mm, and 37.0 mm
<b>Tip Diameters</b>	2.0 mm–3.3 mm tip
<b>Tail Diameters</b>	2.5 mm–4.0 mm tail
<b>Guide Wire</b>	WS-1606DT
<b>Guide Wire Diameter</b>	.062" (1.6 mm)
<b>Hex Driver</b>	HDF-1500 and HD-2000
<b>Hex Size</b>	1.5 mm and 2.0 mm
<b>Drill Part No.</b>	WS-1606DT (optional drill sizes 14.0 mm–22.0 mm)  ATF-024 (14.0 mm, 24.0 mm, and 30.0 mm) ATF-032 (27.0 mm and 32.0 mm) ATF-037 (37.0 mm)

<b>Acutrak® Fusion Hammertoe</b>	<b>Acutrak Fusion Hammertoe</b>  <b>30 mm; 6 mm nose</b>
<b>Acutrak Part No.</b>	30-0033-S
<b>Lengths</b>	30 mm; 6 mm nose
<b>Tip Diameters</b>	1.0 mm tip
<b>Threaded Tip Diameters</b>	2.0 mm threaded tip
<b>Tail Diameters</b>	2.5 mm tail
<b>Guide Wire</b>	WS-1606DT
<b>Guide Wire Diameter</b>	.062" (1.6 mm)
<b>Drill Part No.</b>	ATF-024
<b>Driver Part No.</b>	80-0094

<b>Acutrak®</b>	<b>Acutrak Mini</b>  <b>8.0 mm</b>	<b>Acutrak Standard</b>  <b>12.5 mm</b>	<b>Acutrak 4/5</b>  <b>25.0 mm</b>	<b>Acutrak Plus</b>  <b>35.0 mm</b>	<b>Acutrak 6/7</b>  <b>40.0 mm</b>
<b>Acutrak Part No.</b>	ATM-XXX-S	AT-1XXX-S	AM-00XX-S	AP-00XX-S	AP-67XX-S or AP-671XX-S
<b>Lengths</b>	8.0 mm–26.0 mm (2.0 mm increments)	12.5 mm–30.0 mm (2.5 mm increments)	25.0 mm–50.0 mm (5.0 mm increments)	35.0 mm–80.0 mm (5.0 mm increments)	40.0 mm–120.0 mm (5.0 mm increments)
<b>Tip Diameters</b>	2.8 mm tip	3.3 mm tip	4.0 mm tip	5.2 mm tip	6.0 mm tip
<b>Tail Diameters</b>	3.1 mm–3.6 mm tail	3.8 mm–4.6 mm tail	5.0 mm tail	6.5 mm tail	7.5 mm tail
<b>Guide Wire</b>	WS-0906ST	WS-1106ST	WS-1407ST	WS-1609ST	WS-2408ST
<b>Guide Wire Diameter</b>	.035" (0.9 mm)	.045" (1.1 mm)	.054" (1.4 mm)	.062" (1.6 mm)	.094" (2.4 mm)
<b>Hex Driver</b>	HD-1509	HD-2011	HD-2515 or HP-2515	HD-3016 or HP-3016	HT-4000 (may need T-Handle MS-8043)
<b>Hex Size</b>	1.5 mm	2.0 mm	2.5 mm	3.0 mm	4.0 mm

### Acutrak® Fusion Dimensions

Screws	ATF-140-S	ATF-160-S	ATF-180-S	ATF-200-S	ATF-220-S	ATF-240-S	ATF-270-S	30-0033-S	ATF-320-S	ATF-370-S
<b>Length</b>	14.0 mm	16.0 mm	18.0 mm	20.0 mm	22.0 mm	24.0 mm	27.0 mm	30.0 mm	32.0 mm	37.0 mm
<b>Tip</b>	2.3 mm	2.2 mm	2.2 mm	2.1 mm	2.0 mm	2.0 mm	2.8 mm	2.0 mm	2.8 mm	3.3 mm
<b>Tail</b>	2.5 mm	2.5 mm	2.5 mm	2.5 mm	2.5 mm	2.5 mm	3.5 mm	2.5 mm	3.5 mm	4.0 mm
<b>Tail Shaft Dia.</b>	N/A	N/A	N/A	N/A	N/A	N/A	2.5 mm	N/A	2.5 mm	3.0 mm
<b>Drills</b>	ATF-024	ATF-024	ATF-024	ATF-024	ATF-024	ATF-024	ATF-032	ATF-024	ATF-032	ATF-037
<b>Trocar Tip Lgth.</b>	16.3 mm	16.3 mm	16.3 mm	16.3 mm	16.3 mm	16.3 mm	6.7 mm	16.3 mm	6.7 mm	7.1 mm
<b>Trocar Tip Dia.</b>	1.6 mm	1.6 mm	1.6 mm	1.6 mm	1.6 mm	1.6 mm	1.6 mm	1.6 mm	1.6 mm	1.6 mm
<b>Drill Tip Dia.</b>	2.0 mm	2.0 mm	2.0 mm	2.0 mm	2.0 mm	2.0 mm	2.0 mm	2.0 mm	2.0 mm	2.4 mm
<b>Drill Tail Dia.</b>	2.2 mm	2.2 mm	2.2 mm	2.2 mm	2.2 mm	2.2 mm	3.1 mm	2.2 mm	3.1 mm	3.4 mm
<b>Opt. K-wire for Drill</b>	.062"	.062"	.062"	.062"	.062"	N/A	N/A	N/A	N/A	N/A
<b>1st Band Dia.</b>	26.7 mm	26.7 mm	26.7 mm	26.7 mm	26.7 mm	26.7 mm	30.0 mm	26.7 mm	30.0 mm	35.0 mm
<b>A</b>	29.0 mm	29.0 mm	29.0 mm	29.0 mm	29.0 mm	29.0 mm	32.0 mm	29.0 mm	32.0 mm	38.0 mm
<b>B</b>	32.0 mm	32.0 mm	32.0 mm	32.0 mm	32.0 mm	32.0 mm	35.0 mm	32.0 mm	35.0 mm	42.0 mm
<b>C</b>	35.0 mm	35.0 mm	35.0 mm	35.0 mm	35.0 mm	35.0 mm	37.0 mm	35.0 mm	37.0 mm	45.0 mm
<b>D</b>	38.0 mm	38.0 mm	38.0 mm	38.0 mm	38.0 mm	38.0 mm	41.0 mm	38.0 mm	41.0 mm	50.0 mm

### Drill Tip and Taper End Diameters

	Regular Drill Part No.	Quick Release Part No.	Dense Bone Drill	Tip Diameter	Taper End Diameter
<b>Acutrak Mini</b>	ATM-078 ATM-099	80-0295		2.3 mm 2.3 mm	3.4 mm 3.4 mm
<b>Acutrak Standard</b>	AT-7032 AT-7044	80-0493 80-0296	Dense	2.5 mm 2.8 mm	4.7 mm 4.7 mm
<b>Acutrak 4/5</b>	AM-5010 AM-5014	80-0494 80-0297	Dense	3.0 mm 3.1 mm	4.4 mm 4.6 mm
<b>Acutrak Plus</b>	AP-0100 AP-0104	80-0495 80-0298	Dense	4.0 mm 4.4 mm	5.6 mm 6.0 mm
<b>Acutrak 6/7</b>	AP-67011 AP-67012 AP-67013 AP-67014 AP-67015 AP-67016	80-0496 80-0497 80-0498 80-0477 80-0478 80-0479	Dense Dense Dense	4.4 mm 4.6 mm 4.6 mm 5.1 mm 5.3 mm 5.3 mm	6.7 mm 6.7 mm 6.7 mm 6.7 mm 6.7 mm 6.7 mm

Drill Length from tip to A, B, C, D Marking

